

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

List of contents

| Lot number | Item | 480 wells | 4800 wells |
|------------|--|-----------|------------|
| _____ | Antibody-coated 96-well microtiter plates | 5 | 50 solid |
| _____ | Peroxidase enzyme conjugate, concentrated | 0.550 mL | 1 X 5.5 mL |
| _____ | RUB6, enzyme conjugate diluent | 55 mL | 1 X 550 mL |
| _____ | TMB substrate solution | 60 mL | 550 mL |
| _____ | Positive control | 1 | 5 |
| | The above items should be stored at 2 - 8 °C | | |
| _____ | PBST wash buffer | 7 | 3 X 110 g |
| | The above items should be stored at room temperature (18 - 30 °C). | | |

Materials required, but not provided

Some of the items in the list below may be necessary depending on the type of samples and the method necessary to process the samples. Please refer to sample preparation section for guidance.

- Distilled or purified water
- Paper towels
- Micropipette
- Micropipette tips
- Airtight container for incubations
- Scissors, marker, timer
- Single seed and leaf extraction equipment.
 - Seed press or seed crusher and plate
 - Agdia sample mesh bag (ACC 00930) and rubber mallet
 - Agdia sample mesh bag (ACC 00930) and marker with bag stand
 - Mortar and pestle
 - Micro tube and pestle with tube rack
 - Graduated cylinder
 - Analytical balance
 - Micro tubes and tube rack

Storing the reagents

Store all kit components at the recommended temperature (above) to assure their full shelf life. Each ELISA plate pouch contains a desiccant packet. Keep the plate or unused testwells sealed in the pouch with the desiccant and store in the refrigerator (2 - 8 °C) between uses. Allow the components of the kit to warm to room temperature for about 30 minutes before using.

Safety

Prevent direct skin and eye contact with, or ingestion of, kit components. Obtain medical attention in case of accidental ingestion of kit components. Always wash hands thoroughly after using the kit. It is recommended that gloves be worn when handling the enzyme conjugate solution.

Intended Use

This kit has been validated and approved by Dow AgroSciences® for the detection of the HERCULEX® RW trait. This test system can be used to test individual corn seed and corn leaf.

HERCULEX® RW is a registered trademark of Dow AgroSciences LLC.

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

Test Principle

The test system for Bt-Cry34Ab1 is a direct DAS ELISA. Polyclonal antibodies specific to Bt-Cry34Ab1 are coated to the testwells of a microplate. An enzyme conjugate solution has been included in this kit containing monoclonal antibodies specific to Bt-Cry34Ab1 conjugated to a peroxidase enzyme. Enzyme conjugate is added to the testwells followed by sample extracts. If Bt-Cry34Ab1 is present in the sample, it is bound by the antibodies and captured on the microplate. The plate is then washed to remove any unbound enzyme conjugate and sample. Finally, a substrate is added to the microplate. If peroxidase is present, a color will be produced signifying the presence of Bt-Cry34Ab1. The color reaction can be measured with a plate reader or observed visually.

Please read these instructions carefully before performing the test.

Limitations

The following is a description of factors that could limit test performance or interfere with proper test results.

Buffers: Prepare only the amount of 1X buffers needed for the day. Dilute only the amount of enzyme conjugate necessary at the time of each test run. Do not store 1X buffers.

Samples: This test has been evaluated in corn only.

Sample Extraction Buffer: The Bt-Cry34Ab1 ELISA must be used with 1X PBST wash buffer for optimal results. Do not use sample extraction buffers used with other ELISA kits.

Sample Dilution: ELISA performance is very dependent on the proper sample dilution (tissue weight in g: buffer volume in mL).

Expiration: Test should be used within 1 year from date of purchase.

Storage: Test results may be weak or the test may fail if the storage instructions are not followed properly.

Timing: Please follow as closely as possible the recommended incubation times. Timings for each sample type have been optimized to give the best results for both negative and positive samples. **Note: Please follow tables provided for extraction and incubation times based upon tissue type. Not adhering to these exact times will interfere with achieving proper test results.**

Technical service

If you have any questions about using this kit, please contact Agdia, Inc. Monday – Friday by phone (574-264-2014 or 800-622-4342) or by email (info@agdia.com).

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

Preparing for the test

Familiarize yourself with the kit components and check that all components are present in the kit.

Prepare buffer

Concentrated PBST is used as wash buffer and sample extraction buffer. PBST is supplied as either 20X concentrate or as a powder.

20X concentrate Prepare 1X PBST wash buffer by diluting one 20X pouch of PBST wash buffer with 950 mL of distilled water.

powder Prepare 1X buffer by dissolving PBST buffer powder in distilled water according to the table below:

| | |
|-----------------|--------|
| Buffer powder | 5 g |
| Distilled water | 500 mL |

Prepare controls

Reconstitute lyophilized positive control with 2.0 mL of prepared 1X PBST sample extraction buffer per bottle.

Make control aliquots

After preparing the positive control, divide into aliquots, each sufficient for one use. Dispense aliquots into tubes that can be securely capped. If you will be using a control in one well each time you run the test, prepare 120 μ L aliquots. If you will be using a control in two wells, prepare 220 μ L aliquots. Each aliquot should be sufficient for the tests to be run plus a small additional volume to assure easy dispensing.

Control aliquots must be stored frozen (-10 to -30 °C freezer or household freezer). Do not thaw until just before use. At the time of each test run, remove from storage only the aliquots that will be used. Allow the tubes to thaw and mix the contents thoroughly. At the time you add sample extracts to testwells, add the same volume of control to the appropriate control wells.

Do not refreeze controls.

Prepare testwells

If you will be using less than a full 96-well plate, remove any unused strips and seal them in the foil pouch with the desiccant. Using a permanent marker, number the strips in case a strip becomes separated from the frame.

Prepare a humid box by lining an airtight container with a wet paper towel. Keeping testwells in a humid box during incubation will help prevent samples from evaporating.

Make a copy of the loading diagram and record the locations of your samples and controls. We recommend that you use a buffer well, negative control well and positive control well on each plate each time you run the test.

Preparing Single Leaf and Seed Samples

Leaves, seedlings, or seeds must be ground and diluted in 1X PBST sample extraction buffer. For best results, samples should be diluted in 1X PBST buffer according to the ratios and times listed in the table below.

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

Individual leaves

Sample grinding in Agdia sample mesh bags.



For leaf samples use Agdia's sample mesh bags, a clean mortar and pestle, or any other grinding device that can break up leaf tissues and prevent contamination between samples. This instruction describes the use of Agdia mesh sample bags as the extraction method.

Determine the weight of the leaf sample. Place the leaf sample between the mesh linings of the bag. Rub the bag with a marker to completely crush the sample and to mix the contents uniformly. Make the first dilution by adding 1X PBST buffer to the bag containing the crushed leaf according to a 1:50 leaf to buffer ratio. Massage the bag by hand for a few seconds to ensure good extraction. Let the extract sit for 3 minutes before making the second dilution.

Make a second dilution of the ground sample at a 1:3 extract to buffer ratio by transferring a small volume of the sample extract from the bag and adding it to a tube containing 1X PBST buffer.

| Leaf Tissue | Sample to buffer ratio | Extraction Time |
|--|------------------------|-----------------|
| 1 st dilution | 1:50 g:mL | 3 minutes |
| 2 nd dilution | 1:3 mL:mL | 0 minutes |
| Example | | |
| 1 st dilution 1:50 <ul style="list-style-type: none">• Add 3 mL of 1X PBST to a ground leaf sample weighing 0.06 g.• Mix and allow the extract to sit for 3 minutes. | | |
| 2 nd dilution 1:3 <ul style="list-style-type: none">• Transfer 0.5 mL of extract to a tube containing 1.5 mL of 1X PBST. | | |

Single seeds

Single seeds can be crushed in a seed press, seed crusher or sample mesh bag and rubber mallet. Wash and rinse the grinding equipment between samples.

Determine the weight of the seed. Crush the seed and add the appropriate volume of 1X PBST buffer according to a 1:50 seed to buffer ratio (tissue weight in g: buffer volume in mL). Mix or massage for a few seconds to ensure good extraction. Let the extract sit for 3 minutes before making the next dilution in 1X PBST.

Make the second dilution of the ground sample at a 1:3 extract to buffer ratio by transferring a small volume of the sample extract into a tube of 1X PBST. Mix well before transferring 100 µL of the final sample extract to the ELISA plate.

| Seed Tissue | Sample to buffer ratio | Extraction Time |
|--|------------------------|-----------------|
| 1 st dilution | 1:50 g:mL | 3 minutes |
| 2 nd dilution | 1:3 mL:mL | 0 minutes |
| Example | | |
| 1 st dilution 1:50 <ul style="list-style-type: none">• Add 15 mL of 1X PBST to a ground seed weighing 0.3 g• Mix and allow the extract to sit for 3 minutes. | | |
| 2 nd dilution 1:3 <ul style="list-style-type: none">• Transfer 0.5 mL of extract to a tube containing 1.5 mL of 1X PBST. | | |

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

Test Procedure

- 1. Prepare enzyme conjugate**

The enzyme conjugate is concentrated (100X) and must be diluted with RUB6 enzyme conjugate diluent before use. Prior to use gently shake each vial 10 seconds or vortex for 5 seconds before using. Dilute only what is needed for one day.

Add 110 μ L of concentrated enzyme conjugate to 11 mL of RUB6 diluent, this will be sufficient for 1 plate.

Add 1.1 mL of concentrated enzyme conjugate to 110 mL of RUB6 diluent, this will be sufficient for 10 plates.

Mix the enzyme conjugate thoroughly before adding it to the plate.

Any unused conjugate must be stored in the refrigerator and used within two weeks of diluting.
- 2. Add enzyme conjugate**

Dispense 100 μ L of enzyme conjugate per well.
- 3. Dispense samples, controls, and buffer**

Following your loading diagram, dispense 100 μ L of each prepared sample into sample wells. Dispense 100 μ L of positive control into the positive control wells, 100 μ L of negative control into the negative control wells, and 100 μ L of 1X PBST buffer into the buffer wells.

Mix the contents of the wells by gently swirling the plate on the bench-top.
- 4. Incubate plate**

Set the plate inside the humid box and incubate at room temperature for 60 minutes.
- 5. Warm TMB substrate solution**

About 15 minutes before the end of the above incubation step, measure the required amount of TMB substrate needed. Return the remaining TMB substrate to the refrigerator. Allow measured TMB substrate to warm to room temperature. Caution: TMB substrate is light sensitive, extra precautions are necessary to protect it from light sources when warming to room temperature.

You will need 100 μ L of substrate for each testwell you are using. To estimate the volume needed, measure 1 mL for each 8 well strip used. A full plate will require about 10 mL.
- 6. Wash plate**

When the sample incubation is complete, wash the plate. Use a quick flipping motion to dump the wells into a sink or waste container without mixing the contents.

Fill all the wells completely with 1X PBST, and then quickly empty them again. Repeat 7 times.

After washing, hold the frame upside down and tap firmly on a folded paper towel to remove all droplets of wash buffer.

Note: If using an automatic plate washer, please be sure that the machine is at the appropriate settings for washing flat bottom plates.

Bt-Cry34Ab1 ELISA Kit

Qualitative DAS ELISA for the detection of the Bt-Cry34Ab1 transgenic protein
Catalog number: PSP 04500

7. Add TMB substrate solution Add 100 µL of the TMB substrate solution into each well of the plate.
Let the plate incubate for 30 minutes. Keep testwells away from strong light.
8. Evaluate results Measure O.D.'s on a plate reader at 650 nm. Air bubbles which are present at the time of reading can alter results, if in the light path. Agdia recommends that bubbles be eliminated prior to reading.
- Wells in which a blue color develops indicate positive results. Wells in which there is no significant color development indicate negative results. Test results are valid only if positive control wells give a positive result and buffer wells remain colorless.
- If either control well does not show the appropriate color, please repeat the test procedure. If the problem persists, contact Agdia for further assistance.
- Optional: Stop solution may be added to terminate the peroxidase/TMB reaction. Add 1N HCl at 100 µL per testwell. Read at 450 nm up to one hour after addition.

Buffer Formulations

The following buffer is a standard part of your kit. This formulation is for reference only.

PBST Buffer (Wash Buffer) (1X)

Dissolve in distilled water to 1000 mL:

| | |
|--|--------|
| Sodium chloride | 8.0 g |
| Sodium phosphate, dibasic (anhydrous) | 1.15 g |
| Potassium phosphate, monobasic (anhydrous) | 0.2 g |
| Potassium chloride | 0.2 g |
| Tween-20 | 0.5 g |

Adjust pH to 7.4

HERCULEX[®] RW is a registered trademark of Dow AgroSciences LLC.

Date _____ Test _____

Test performed by _____

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | | | | | | | | | | | | |
| B | | | | | | | | | | | | |
| C | | | | | | | | | | | | |
| D | | | | | | | | | | | | |
| E | | | | | | | | | | | | |
| F | | | | | | | | | | | | |
| G | | | | | | | | | | | | |
| H | | | | | | | | | | | | |

