MedChemExpress

®

'N´

CI

Q

юн

CI

Bendamustine

Cat. No.:	HY-13567
CAS No.:	16506-27-7
Molecular Formula:	C ₁₆ H ₂₁ Cl ₂ N ₃ O ₂
Molecular Weight:	358.26
Target:	DNA Alkylator/Crosslinker; Apoptosis
Pathway:	Cell Cycle/DNA Damage; Apoptosis
Storage:	-20°C, protect from light
	* In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)

SOLVENT & SOLUBILITY

In Vitro	DMSO : 100 mg/mL (279.13 mM; Need ultrasonic)						
	Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg		
		1 mM	2.7913 mL	13.9563 mL	27.9127 mL		
		5 mM	0.5583 mL	2.7913 mL	5.5825 mL		
		10 mM	0.2791 mL	1.3956 mL	2.7913 mL		
	Please refer to the so	lubility information to select the ap	opropriate solvent.				
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (6.98 mM); Clear solution						
	2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (6.98 mM); Clear solution						
	3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (6.98 mM); Clear solution						

BIOLOGICAL ACTIVITY				
BIOLOGICAL ACTIVITY				
Description	Bendamustine (SDX-105 free base), a purine analogue, is a DNA cross-linking agent. Bendamustine activates DNA-damage stress response and apoptosis. Bendamustine has potent alkylating, anticancer and antimetabolite properties ^[1] .			
IC ₅₀ & Target	DNA Alkylator/Crosslinker ^[1]			
In Vitro	Bendamustine is a DNA cross-linking agent that causes DNA breaks, with alkylating and antimetabolite properties. Bendamustine uniquely regulates apoptosis pathways and DNA repair pathways in non-Hodgkin's lymphoma cells. Bendamustine (50 μM) induces p21 (Cip1/Waf1) and NOXA genes, and increases the expression of p53 in SU-DHL-1 cells. Bendamustine (25 μM) blocks mitotic checkpoints and cuases mitotic catastrophe ^[1] .			

	Bendamustine reduces the viability of multiple myeloma (MM) cell lines, such as RPMI-8226 and 8226-LR5 cells, with IC ₂₅ s of 101.8 μM and 585.5 μM after 24 h incubation, and 51.7 and 374.3 μM after 48 h incubation, respectively. Bendamustine induces a specific caspase-dependent MM cell death and inhibits the spindle-assembly checkpoint ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
In Vivo	Bendamustine (25 mg/kg, i.v.) shows potent inhibition on the growth of tumor cells by 91%, 99% and 95% for DoHH-2, Granta 519 and RAMOS models, respectively. Moreover, the antitumor effect of Bendamustine is enhanced by rituximab in DoHH-2 and RAMOS models, but not in Granta 519 model ^[3] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Cell Assay^[2]

Cytotoxicity of both Bendamustine and melphalan on multiple myeloma (MM) cells is calculated as inhibition of cell viability by measuring the percentage of cell survival by MTS assay. Briefly, cells (1×10^4 /well) are seeded in 96-well plates with increasing concentrations of the drug and analyzed after 24, 48, 72 and 96 h of incubation. To this end, 1 µg/mL of MTS solution is added to each well and, after 1 h at 37 °C, the dark blue formazan crystals are dissolved by isopropanol 1 N and HCl (24:1, vol/vol). Finally, the absorbance is measured at 490 nm in a 96-well plate reader. Cell survival is estimated as the percentage of the absorbance of untreated controls and each test is performed in triplicate. The inhibitory concentrations 50 (IC₅₀) and 25 (IC₂₅) of each drug, being the amount able to reduce cell growth to 50% and 25%, respectively, of that of untreated control cells, are calculated, and the tests are performed in parallel using equitoxic concentrations of Bendamustine and melphalan. The relative resistance index (RRI) is expressed as the ratio of the IC₅₀ of 8226-LR5 to the IC₅₀ of RPMI-8226 cells^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration ^[3]

C.B.-17 scid mice (DoHH-2, Granta 519) or C.B.-17 scid-bg mice (SuDHL-4, RAMOS) are inoculated with 1×10^{6} (DoHH-2, RAMOS), 3×10^{6} (SuDHL-4) or 5×10^{6} (Granta 519) cells s.c. in the right flank. For flank xenografts, inoculation volume is 0.2 mL consisting of a 50:50 mixture of cells in growth medium and Matrigel. Tumour volume is estimated by two to three weekly measurements of the length and width of the tumour by electronic calipers and applying the following equation: $V=L\times W^{2}/2$. Tumours are allowed to reach approximately 250 mm³, and mice are size-matched (day 0) into treatment and control groups. For systemic Granta 519 tumour models, 2×10^{6} cells are injected via the tail vein in 0.1 mL volume of cell medium on day 0, and treatment is initiated on day 14. All animals are ear-tagged and monitored individually throughout the experiment. Navitoclax is administered by oral gavage once daily in a mixture of Phosal 50PG : PEG400 : ethanol. Bendamustine and rituximab are administered i.v. at 25 mg/kg and 10 mg/kg, respectively, on day 1. Navitoclax is administered i.v. at 25 mg/kg and 10 mg/kg, respectively, on day 1. Navitoclax is include loss of ambulation, laboured breathing or weight loss > 20% mean body weight per cage^[3]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- J Mol Med (Berl). 2019 Aug;97(8):1183-1193.
- J Biomed Res. 2017 0(0): 1-12.
- Patent. US20160222465A1.

See more customer validations on www.MedChemExpress.com

Mice^[3]

REFERENCES

[1]. Leoni LM, et al. Bendamustine (Treanda) displays a distinct pattern of cytotoxicity and unique mechanistic features compared with other alkylating agents. Clin Cancer Res. 2008 Jan 1;14(1):309-17.

[2]. Cives M, et al. Bendamustine overcomes resistance to melphalan in myeloma cell lines by inducing cell death through mitotic catastrophe. Cell Signal. 2013 May;25(5):1108-17.

[3]. Ackler S, et al. Navitoclax (ABT-263) and bendamustine ± rituximab induce enhanced killing of non-Hodgkin's lymphoma tumours in vivo. Br J Pharmacol. 2012 Oct;167(4):881-91.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA