DOTAP chloride

MedChemExpress

®

HY-112754A			
132172-61-3			
C ₄₂ H ₈₀ CINO ₄			
698.54			
Liposome			
Metabolic Enzyme/Protease			
-20°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)			
	HY-112754A 132172-61-3 C ₄₂ H ₈₀ ClNO ₄ 698.54 Liposome Metabolic Enzyme/Protease -20°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)		

SOLVENT & SOLUBILITY

In Vitro	Ethanol : 130 mg/mL (186.10 mM; Need ultrasonic) DMSO : 25 mg/mL (35.79 mM; ultrasonic and warming and heat to 60°C)					
	Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg	
		1 mM	1.4316 mL	7.1578 mL	14.3156 mL	
		5 mM	0.2863 mL	1.4316 mL	2.8631 mL	
		10 mM	0.1432 mL	0.7158 mL	1.4316 mL	
	Please refer to the so	lubility information to select the app	propriate solvent.			
In Vivo	1. Add each solvent Solubility: ≥ 2.5 m	one by one: 10% DMSO >> 40% PEC g/mL (3.58 mM); Clear solution	G300 >> 5% Tween-80) >> 45% saline		
	2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: 2.5 mg/mL (3.58 mM); Suspended solution; Need ultrasonic					
	3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (3.58 mM); Clear solution					

Description	DOTAP chloride is a useful and effective cationic lipid for transient and stable transfection DNA (plasmids, bacmids) and modified nucleic acids (antisense oligonucleotides) with out the use of helper lipid ^[1] .				
IC ₅₀ & Target	IC50: Cationic Lipid ^[1]				
In Vitro	liposomes are reorganized by associating with DNA. The transfection efficiency of DOTAP liposomes is mainly influenced by lipid composition and cell type, but not by size or zeta potential ^[1] .				

 $\sim\sim$

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Adv Healthc Mater. 2022: 2202412.
- Pharmaceutics. 2023 Apr 30, 15(5), 1379.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Kim BK, et al. DOTAP/DOPE ratio and cell type determine transfection efficiency with DOTAP-liposomes. Biochim Biophys Acta. 2015 Oct;1848(10 Pt A):1996-2001.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA