# **Screening Libraries**



## Gilteritinib

Cat. No.: HY-12432 1254053-43-4 CAS No.: Molecular Formula:  $C_{29}H_{44}N_8O_3$ Molecular Weight: 552.71

Target: FLT3; TAM Receptor

Pathway: Protein Tyrosine Kinase/RTK Storage: 4°C, stored under nitrogen

\* In solvent: -80°C, 6 months; -20°C, 1 month (stored under nitrogen)

**Product** Data Sheet

### **SOLVENT & SOLUBILITY**

In Vitro

Ethanol: 100 mg/mL (180.93 mM; ultrasonic and adjust pH to 2 with HCl)

DMSO: 2 mg/mL (3.62 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	1.8093 mL	9.0463 mL	18.0927 mL
	5 mM	0.3619 mL	1.8093 mL	3.6185 mL
	10 mM	0.1809 mL	0.9046 mL	1.8093 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

1. Add each solvent one by one: 50% PEG300 >> 50% saline Solubility: 10 mg/mL (18.09 mM); Suspended solution; Need ultrasonic

### **BIOLOGICAL ACTIVITY**

Description	$ \textit{Gilteritinib (ASP2215)} \ is \ a \ potent \ and \ ATP-competitive \ FLT3/AXL \ inhibitor \ with \ IC_{50}s \ of \ 0.29 \ nM/0.73 \ nM, \ respectively. $
IC <sub>50</sub> & Target	IC50: 0.29 nM (FLT3) <sup>[1]</sup> IC50: 0.35 nM (LTK), 0.73 nM (AXL), 1.2 nM (EML4-ALK), 230 nM (c-KIT) <sup>[2]</sup>
In Vitro	Of the 78 tyrosine kinases tested, Gilteritinib (ASP2215) inhibits FLT3, leukocyte tyrosine kinase (LTK), anaplastic lymphoma kinase (ALK), and AXL kinases by over 50% at 1 nM with an IC <sub>50</sub> value of 0.29 nM for FLT3, approximately 800-fold more potent than for c-KIT <sup>[1]</sup> . Gilteritinib inhibits the activity of eight of the 78 tested kinases by over 50% at concentrations of either 1 nM (FLT3, LTK, ALK, and AXL) or 5 nM (TRKA, ROS, RET, and MER). The IC <sub>50</sub> s are 0.29 nM for FLT3 and 0.73 nM for AXL. Gilteritinib inhibits FLT3 at an IC <sub>50</sub> that is approximately 800-fold more potent than the concentration required to inhibit c-

KIT (230 nM). The antiproliferative activity of Gilteritinib is evaluated against MV4-11 and MOLM-13 cells, which

endogenously express FLT3-ITD. After 5 days of treatment, Gilteritinib inhibits the growth of MV4-11 and MOLM-13 cells with mean IC<sub>50</sub>s of 0.92 nM (95% CI: 0.23-3.6 nM) and 2.9 nM (95% CI: 1.4-5.8 nM), respectively. Growth suppression of MV4-11

Page 1 of 3

cells is accompanied by inhibition of FLT3 phosphorylation. Relative to vehicle control cells, phosphorylated FLT3 levels are 57%, 8%, and 1% after 2 h of treatment with 0.1 nM, 1 nM, and 10 nM Gilteritinib, respectively. In addition, doses as low as 0.1 nM or 1 nM result in the suppression of phosphorylated ERK, STAT5, and AKT, all of which are downstream targets of FLT3 activation. To investigate the effects of Gilteritinib on AXL inhibition, MV4-11 cells that expressed exogenous AXL are treated with Gilteritinib. At concentrations of 1 nM, 10 nM, and 100 nM for 4 h, Gilteritinib treatment decreases phosphorylated AXL levels by 38%, 29%, and 22%, respectively<sup>[2]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### In Vivo

In MV4-11 xenografted-mice, the concentration of Gilteritinib (ASP2215) in tumors is more than 20-fold higher than that in plasma with oral administration of Gilteritinib at 10 mg/kg for 4 days. Treatment of Gilteritinib for 28 days results in dose-dependent inhibition of MV4-11 tumor growth and induces complete tumor regression at more than 6 mg/kg. Further, Gilteritinib decreases tumor burden in bone marrow and prolonged the survival of mice intravenously transplanted with MV4-11 cells<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### **PROTOCOL**

### Kinase Assay [2]

The kinase inhibitory activity of Gilteritinib is tested against a panel of 78 tested kinases using ATP concentrations that are approximately equal to the  $K_m$  value for each kinase in a TK-ELISA or off-chip mobility shift assay. Initially, two concentrations of Gilteritinib (1 nM and 5 nM) are tested to assess each compound's inhibitory effect on TK activity. Further studies are then conducted using a dose range of Gilteritinib to determine  $IC_{50}$  values for kinases in which activity is inhibited by >50% with 1 nM Gilteritinib as well as for c-KIT. TK-ELISA and MSA assays are used to conduct  $IC_{50}$  studies for FLT3, LTK, AXL, and c-KIT; the HTRF KinEASE-TK assay is performed to assess the  $IC_{50}$  value of echinoderm microtubule-associated protein-like 4-ALK (EML4-ALK)<sup>[2]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### Cell Assay [2]

The effect of Gilteritinib on MV4-11 and MOLM-13 cells is assessed using the CellTiter-Glo Luminescent Cell Viability Assay. Subsequent studies are conducted to examine the effect of Gilteritinib and Quizartinib on Ba/F3 cells expressing either FLT3-ITD, FLT3-D835Y, FLT3-ITD-D835Y, FLT3-ITD-F691 L, or FLT3-ITD-F691I. MV4-11 and MOLM-13 cells are treated with DMSO or increasing concentrations of Gilteritinib (0.01, 0.1, 1, 10, and 100 nM) for 5 days, and cell viability is measured using CellTiter-Glo<sup>[2]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

# Animal Administration [1]

### Mice<sup>[1]</sup>

Antitumor activity is evaluated in nude mice transplanted with MV4-11 AML cells. The pharmacokinetics in xenografted mice is also investigated. MV4-11 xenografted-mice are treated with oral administration of Gilteritinib at 10 mg/kg for 4 days. Treatment of Gilteritinib for 28 days results in dose-dependent inhibition of MV4-11 tumor growth and induces complete tumor regression at more than 6 mg/kg $^{[1]}$ .

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### **CUSTOMER VALIDATION**

- Science. 2017 Dec 1;358(6367):eaan4368.
- Cancer Discov. 2023 Apr 3;CD-22-0411.
- Sci Adv. 2022 Sep 16;8(37):eabp9005.
- Blood Cancer J. 2022 Jan 11;12(1):5.
- Haematologica. 2018 Nov;103(11):1862-1872.

See more customer validations on <u>www.MedChemExpress.com</u>						
REFERENCES						
[1]. ASP2215, a novel FLT3/AXL inhibitor: Preclinical evaluation in acute myeloid leukemia (AML). 2014 ASCO Annual Meeting.						
[2]. Mori M, et al. Gilteritinib, a FLT3/AXL inhibitor, shows antileukemic activity in mouse models of FLT3 mutated acute myeloid leukemia. Invest New Drugs. 2017 Oct;35(5):556-565.						
	_					
Caution: Product has not been fully validated for medical applications. For research use only.						
Tel: 609-228-6898 Fax: 609-228-5909 E-mail: tech@MedChemExpress.com  Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA						