


## Reldesemtiv

<b>Cat. No.:</b>	HY-109121		
<b>CAS No.:</b>	1345410-31-2		
<b>Molecular Formula:</b>	C <sub>19</sub> H <sub>18</sub> F <sub>2</sub> N <sub>6</sub> O		
<b>Molecular Weight:</b>	384.38		
<b>Target:</b>	Others		
<b>Pathway:</b>	Others		
<b>Storage:</b>	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	2 years
		-20°C	1 year



### SOLVENT & SOLUBILITY

<b>In Vitro</b>	DMSO : 250 mg/mL (650.40 mM; Need ultrasonic)				
		Solvent Concentration	Mass 1 mg	5 mg	10 mg
	<b>Preparing Stock Solutions</b>	1 mM	2.6016 mL	13.0080 mL	26.0159 mL
		5 mM	0.5203 mL	2.6016 mL	5.2032 mL
10 mM		0.2602 mL	1.3008 mL	2.6016 mL	
Please refer to the solubility information to select the appropriate solvent.					
<b>In Vivo</b>	<ol style="list-style-type: none"> <li>Add each solvent one by one: 10% DMSO &gt;&gt; 40% PEG300 &gt;&gt; 5% Tween-80 &gt;&gt; 45% saline Solubility: ≥ 2.08 mg/mL (5.41 mM); Clear solution</li> <li>Add each solvent one by one: 10% DMSO &gt;&gt; 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.08 mg/mL (5.41 mM); Clear solution</li> <li>Add each solvent one by one: 10% DMSO &gt;&gt; 90% corn oil Solubility: ≥ 2.08 mg/mL (5.41 mM); Clear solution</li> </ol>				

### BIOLOGICAL ACTIVITY

<b>Description</b>	Reldesemtiv (CK-2127107) is a selective, orally active and next-generation fast skeletal muscle troponin activator (FSTA). Reldesemtiv selectively activates fast skeletal myofibrils with an EC <sub>50</sub> of 3.4 μM. Reldesemtiv increases exercise performance in a heart failure model <sup>[1]</sup> .
<b>In Vitro</b>	Reldesemtiv selectively binds to and sensitizes the fast skeletal troponin complex to calcium <sup>[1]</sup> . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## In Vivo

Reldesemtiv (10 mg/kg; p.o.) significantly improves rotarod performance in exercise-intolerant LAD-HF rats<sup>[1]</sup>.

Reldesemtiv (10 mg/kg; i.v.) increases Ca<sup>2+</sup> sensitivity in fast skeletal muscle fibers and increases muscle force in response to submaximal stimulation frequencies in LAD-HF rats<sup>[1]</sup>.

Reldesemtiv increases sarcomeric Ca<sup>2+</sup> sensitivity in skinned LAD-HF skeletal fibers<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Model:	Female Sprague-Dawley rats (Sham and left anterior descending coronary artery (LAD) heart failure (LAD-HF)) <sup>[1]</sup>
Dosage:	10 mg/kg
Administration:	p.o.
Result:	Increased performance time by approximately 150% in LAD-HF rats (116±22 versus 283±47 seconds).

## REFERENCES

[1]. Hwee DT, et al. The small-molecule fast skeletal troponin activator, CK-2127107, improves exercise tolerance in a rat model of heart failure. J Pharmacol Exp Ther. 2015 Apr;353(1):159-68.

[2]. Andrews JA, et al. CK-2127107 amplifies skeletal muscle response to nerve activation in humans. Muscle Nerve. 2018 May;57(5):729-734.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA