TLR7 agonist 2

MedChemExpress

Cat. No.:	HY-103039		
CAS No.:	1642857-69	-9	
Molecular Formula:	$C_{17}H_{16}N_6O_2$		
Molecular Weight:	336.35		
Target:	Toll-like Receptor (TLR)		
Pathway:	Immunology/Inflammation		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	2 years
		-20°C	1 year

R

SOLVENT & SOLUBILITY

		Solvent Mass Concentration	1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM	2.9731 mL	14.8655 mL	29.7309 mL	
	5 mM	0.5946 mL	2.9731 mL	5.9462 mL	
		10 mM	0.2973 mL	1.4865 mL	2.9731 mL

Description	TLR7 agonist 2 is a potent and selective Toll-like Receptor 7 (TLR7) agonist with a LEC of 0.4 $\mu M.$			
IC ₅₀ & Target	LEC: 0.4 μM (TLR7) ^[1]			
In Vitro	TLR7 agonist 2 is a potent and selective Toll-like Receptor 7 (TLR7) agonist with a lowest effective concentration (LEC) of 0.4 μM in HEK293 cell. TLR7 agonist 2 is found to be selective for TLR7 over TLR8 with LEC of >100 μM for human TLR8. TLR7 agonist 2 demonstrates low inhibition across five CYP450 isozymes (IC ₅₀ >10 μM) and is also not a time dependent inhibitor of CYP450 3A4. TLR7 agonist 2 has limited inhibition of the hERG potassium ion channel ³ H-dofetilide binding in vitro (IC ₅₀ >50 μM) ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.			
In Vivo	TLR7 agonist 2 is found to be rapidly cleared in conjunction with our target profile. Both C _{max} and AUC increase less than dose proportionally between 0.3 and 3 mg/kg and more than dose-proportionally between 3 and 10 mg/kg. TLR7 agonist 2 can induce an antiviral interferon stimulated gene (ISG) response without inducing an IFNα response at a low dose. TLR7 agonist 2 also induces a 2.7 log decrease in serum HBV viral load from 0.3 mg/kg, and a maximum 3.1 log decrease is			

 H_2N

li N

Ν

	observed for doses between 1 and 5 mg/kg ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
ROTOCOL	
ell Assay [1]	The ability of TLR7 agonist 2 to activate human TLR7 and/or TLR8 is assessed by using HEK293 cells. Briefly, HEK293 cells are grown in culture medium (DMEM supplemented with 10% FCS and 2 mM Glutamine). Transfected cells are then detached with Trypsin-EDTA, washed in PBS and resuspended in medium to a density of 1.67×10 ⁵ cells/mL. Thirty microliters of cells are then dispensed into each well in 384-well plates, where 10 μL of TLR7-agonist-1 in 4% DMSO is already present. Following 6 hours incubation at 37°C, 5% CO ₂ , the luciferase activity is determined by adding 15 μL of Steady Lite Plus substrate to each well and readout performed on a microplate imager. Lowest effective concentrations (LEC) values are determined for TLR7-agonist-1 ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
nimal	A may so in vivo model is used to domenstrate the initial proof of concent to induce endegenous JENIX. Single eral

AnimalA mouse in vivo model is used to demonstrate the initial proof of concept to induce endogenous IFNα. Single oralAdministration ^[1]administration of 0.3, 1, 3, and 10 mg/kg doses of TLR7 agonist 2 is given to healthy, female, fasted C57Bl/6 mice.
Concentrations of TLR7 agonist 2 and mouse-IFN via ELISA are measured from the plasma and compare to vehicle^[1].
MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

• Research Square Print. 2023 Feb 1.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. McGowan DC, et al. Identification and Optimization of Pyrrolo[3,2-d]pyrimidine Toll-like Receptor 7 (TLR7) Selective Agonists for the Treatment of Hepatitis B. J Med Chem. 2017 Jul 27;60(14):6137-6151.

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech

9 E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA