Proteins

Inhibitors

THZ1 Hydrochloride

Cat. No.: HY-80013A

Molecular Formula: $C_{31}H_{29}Cl_2N_7O_2$ Molecular Weight: 602.51

Target:

Pathway: Cell Cycle/DNA Damage

CDK

Storage: 4°C, sealed storage, away from moisture

* In solvent: -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)

Product Data Sheet

SOLVENT & SOLUBILITY

In Vitro DMSO: 22.5 mg/mL (37.34 mM; Need ultrasonic and warming)

H₂O: < 0.1 mg/mL (ultrasonic; warming; heat to 60°C) (insoluble)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	1.6597 mL	8.2986 mL	16.5972 mL
	5 mM	0.3319 mL	1.6597 mL	3.3194 mL
	10 mM	0.1660 mL	0.8299 mL	1.6597 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.17 mg/mL (3.60 mM); Clear solution
- 2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: 2.17 mg/mL (3.60 mM); Suspended solution; Need ultrasonic
- 3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.17 mg/mL (3.60 mM); Clear solution

BIOLOGICAL ACTIVITY

Description	THZ1 Hydrochloride is a selective and potent covalent CDK7 inhibitor with an IC ₅₀ of 3.2 nM. THZ1 Hydrochloride also inhibits the closely related kinases CDK12 and CDK13 and downregulates MYC expression ^{[1][2]} .			
IC ₅₀ & Target	CDK7	CDK12	CDK13	

3.2 nM (IC₅₀)

In Vitro THZ1 inhibits Jurkat cell and Loucy cell with IC $_{50}$ of 50 nM, and 0.55 nM, respectively. THZ1 demonstrates time-dependent inhibition of CDK7 in vitro and covalent binding of intracellular CDK7. THZ1 (9, 27, 83, 250, 750, and 2500 nM) inhibits CDK12 but at higher concentrations compared to CDK7. THZ1 (1 μ M) irreversibly inhibits RNAPII CTD and CAK phosphorylation. THZ1 (2.5 μ M) irreversibly inhibits RNAPII CTD phosphorylation by covalently targeting a unique cysteine located outside the kinase domain of CDK7 in Hela S3 cells. THZ1 (250 nM) causes decreased cellular proliferation and an increase in apoptotic index with concomitant reduction in anti-apoptotic proteins, most notably MCL-1 and XIAP in T-ALL cell lines^[1]. ?All genotypically-distinct human (hSCLC) cell lines exhibit high sensitivity to THZ1, with an IC₅₀ in the range of 5-20 nM^[3]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

In Vivo

THZ1 (10 mg/kg) demonstrates potent killing of primary chronic lymphocytic leukemia (CLL) cells and anti-proliferative activity against primary TALL cells and in vivo against a human T-ALL xenograft^[1].

?THZ1 (10 mg/kg, i.v.) inhibits tumor growth in a mouse model of human MYCN-amplified NB and shows no toxicity^[4]. ?THZ1 (10 mg/kg, i.p.) completely suppresses oesophageal squamous cell carcinoma tumour growth in vivo without loss of body weight or other common toxic effects^[5].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Kinase Assay [1]

For kinase assays following immunoprecipitation of FLAG-CDK7 protein from HCT116 or FLAG-CDK12 from 293A cellular lysates, cells are first treated with THZ1, THZ1-R, or DMSO for 4 hrs at 37°C. Cells are then harvested by lysis in 50 mM Tris HCl pH 8.0, 150 mM NaCl, 1% NP-40, 5 mM EDTA, and protease/phosphatase cocktails. Exogenous CDK7 or CDK12 proteins are immunoprecipitated from cellular lysates using FLAG antibody- conjugated agarose beads. Precipitated proteins are washed with lysis buffer 6 times, followed by 2 washes with kinase buffer (40 mM Hepes pH 7.5, 150 mM NaCl, 10 mM MgCl₂, 5% glycerol) and subjected to in vitro kinase assays at 30°C for 45 minutes using 1 μ g of the large subunit of RNAPII (RPB1) as substrate and 25 μ M ATP and 10 μ Ci of 32P ATP. Note: no additional inhibitors are added to the kinase reaction mixture, therefore any inhibition results from the activity of the compounds that are added directly to cells (ie-intracellular inhibition). This suggests that the inhibitors are either covalently tethered to the kinase or have strong non-covalent character. Kinase assays using recombinant CDK7/TFIIH/MAT1 are conducted in the manner as described above using 25 ng of CAK complex per reaction. For kinase assays designed to test time-dependent inactivation of CDK7 kinase activity, CAK complex is pre-incubated with indicated concentrations of THZ1, THZ1-R, or DMSO in kinase buffer without ATP for 4 hrs at 30°C prior to being subjected to kinase assay conditions [1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Cell Assay [1]

Jurkat, Loucy, KOPTK1, and DND-41 cell lines are seeded in 384-well microplates at 15% confluency in medium with 5% FBS and penicillin/streptavidin. Cells are treated with THZ1 (2, 10, 50, 250, 1250, and 6250 nM) or DMSO for 72 hrs and cell viability is determined using resazurin^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration [1]

Mice^[1]

Thirty-two NOD-SCIDIL2Rc γ null (NSG) 9-week old female mice are divided into treatment groups based on mean BLI as follows: THZ1 10 mg/kg qD, THZ1 10 mg/kg BID, and vehicle (10% DMSO in D5W) BID (n=10 for all groups). Two mice are excluded, one with the highest and one with the lowest BLI. All treatments are administered via IV injection in the lateral tail vein in a volume of 3.3 μ L/g (non-blinded). Mice are imaged and weighed every 3-5 days. Mice are treated for four weeks and on the final day mice are imaged, dosed and sacrificed approximately 5-6 hrs post dose. Upon sacrifice, blood is collected via cardiac puncture in EDTA tubes; a portion (300 uL) is processed for plasma. Liver and spleen tissues are collected from each mouse with half of each sample flash frozen and half of each sample fixed. Blood plasma and liver samples are processed for pharmacokinetics analysis of THZ1. Spleen tissues are homogenized and lysed and processed for pharmacodynamics analysis of THZ1 target engagement.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Nat Med. 2019 Feb;25(2):292-300.
- Cell. 2018 Sep 20;175(1):171-185.e25.
- Cell. 2017 Sep 7;170(6):1209-1223.e20.
- Science. 2021 Apr 30;372(6541):eaba8490.
- Cancer Discov. 2019 Nov;9(11):1538-1555.

See more customer validations on www.MedChemExpress.com

REFERENCES

- [1]. Kwiatkowski N, et al. Targeting transcription regulation in cancer with a covalent CDK7 inhibitor. Nature. 2014 Jul 31;511(7511):616-20.
- [2]. Jiang YY, et al. Targeting super-enhancer-associated oncogenes in oesophageal squamous cell carcinoma. Gut. 2016 May 10. pii: gutjnl-2016-311818.
- [3]. Christensen CL, et al. Targeting transcriptional addictions in small cell lung cancer with a covalent CDK7 inhibitor. Cancer Cell. 2014 Dec 8;26(6):909-22.
- [4]. Chipumuro, et al. CDK7 inhibition suppresses super-enhancer-linked oncogenic transcription in MYCN-driven cancer. Cell. 2014 Nov 20;159(5):1126-39.?

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA