

ATTO 680 cadaverin

Cat. No.:	HY-D2070
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.

ATTO 680 cadaverin

BIOLOGICAL ACTIVITY

Description

ATTO 680 cadaverin is a cadaverine derivative of ATTO 680, used for carboxylic acid, aldehyde and ketone reactions, maximum excitation/emission wavelength: 681/698 nm.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA