**Proteins** 

# **Screening Libraries**

# **Product** Data Sheet

### A-887826

Cat. No.: HY-100080 CAS No.: 1266212-81-0 Molecular Formula:  $C_{26}H_{29}CIN_4O_3$ Molecular Weight: 480.99

Target: Sodium Channel

Pathway: Membrane Transporter/Ion Channel

-20°C Storage: Powder 3 years

4°C 2 years

-80°C In solvent 6 months

> -20°C 1 month

### **SOLVENT & SOLUBILITY**

In Vitro

DMSO: 250 mg/mL (519.76 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	2.0790 mL	10.3952 mL	20.7905 mL
	5 mM	0.4158 mL	2.0790 mL	4.1581 mL
	10 mM	0.2079 mL	1.0395 mL	2.0790 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.08 mg/mL (4.32 mM); Clear solution

## **BIOLOGICAL ACTIVITY**

Description A-887826 is a potent, selective, oral bioavailable and voltage-dependent Na(v)1.8 sodium channel blocker with an IC50 of 11 nM . A-887826 attenuates neuropathic tactile allodynia in vivo<sup>[1]</sup>.

IC<sub>50</sub> & Target Na<sub>v</sub>1.8

1.1 nM (IC<sub>50</sub>)

In Vitro

A-887826 is approximately 3 fold less potent to block Na(v)1.2, approximately 10 fold less potent to block tetrodotoxinsensitive sodium (TTX-S Na(+)) currents and is >30 fold less potent to block Na(v)1.5 channels<sup>[1]</sup>.

A-887826 potently blocks tetrodotoxin-resistant sodium (TTX-R Na(+)) currents (IC<sub>50</sub>=8 nM) from small diameter rat dorsal root ganglion (DRG) neurons in a voltage-dependent fashion<sup>[1]</sup>.

A-887826 effectively suppresses evoked action potential firing when DRG neurons are held at depolarized potentials and reversibly suppresses spontaneous firing in small diameter DRG neurons from complete Freund's adjuvant inflamed rats<sup>[1]</sup>.

	affecting voltage-depe	A-887826 (100 nM) shifts the mid-point of voltage-dependent inactivation of TTX-R currents by approximately 4mV without affecting voltage-dependent activation and do not exhibit frequency-dependent inhibition <sup>[1]</sup> .  MCE has not independently confirmed the accuracy of these methods. They are for reference only.		
In Vivo	[1].	A-887826 (3-30 mg/kg; p.o.; 1 hour pre-treatment) significantly attenuates tactile allodynia in a rat neuropathic pain model [1].  MCE has not independently confirmed the accuracy of these methods. They are for reference only.		
	Animal Model:	Male Spraguee Dawley rats (200-300 g) <sup>[1]</sup>		
	Dosage:	3 mg/kg, 10 mg/kg, 30 mg/kg		
	Administration:	Oral administration; 1 hour before testing; fourteen days after L5-L6 spinal nerve ligation		
	Result:	Produced a dose-dependent attenuation of tactile allodynia in this experimental pain model.		

### **REFERENCES**

[1]. Zhang XF, et al. A-887826 is a structurally novel, potent and voltage-dependent Na(v)1.8 sodium channel blocker that attenuates neuropathic tactile allodynia in rats. Neuropharmacology. 2010 Sep;59(3):201-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

 $\hbox{E-mail: } tech @ Med Chem Express.com$ 

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA