# R(+)-Methylindazone

Cat. No.: HY-12693 CAS No.: 54197-31-8 Molecular Formula: C<sub>17</sub>H<sub>18</sub>Cl<sub>2</sub>O<sub>4</sub> Molecular Weight: 357.23

Chloride Channel Target:

Pathway: Membrane Transporter/Ion Channel

-20°C Storage: Powder

2 years

3 years

-80°C In solvent 2 years

> -20°C 1 year

**Product** Data Sheet

## **SOLVENT & SOLUBILITY**

In Vitro

DMSO : ≥ 100 mg/mL (279.93 mM)

H<sub>2</sub>O: < 0.1 mg/mL (insoluble)

\* "≥" means soluble, but saturation unknown.

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	2.7993 mL	13.9966 mL	27.9932 mL
	5 mM	0.5599 mL	2.7993 mL	5.5986 mL
	10 mM	0.2799 mL	1.3997 mL	2.7993 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (7.00 mM); Clear solution
- 2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (7.00 mM); Clear solution
- 3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (7.00 mM); Clear solution

### **BIOLOGICAL ACTIVITY**

Description

R(+)-Methylindazone (R(+)-IAA-94) is a potent indanyloxyacetic acid blocker of epithelial chloride channels. R(+)- $Methyl indazone in hibits Nef-sdAb19 \ (single-domain antibody) interaction and binds to negative factor \ (Nef)^{[1][2][3]}.$ 

In Vitro

IAA-94 has been employed in modulating chloride channel function to probe the dynamics and function of the channels. The high affinity of IAA-94 for the chloride channel has been exploited for isolation and reconstitution of these proteins<sup>[1][2]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

# **CUSTOMER VALIDATION**

- Mol Metab. 2023 Aug 19;101794.
- Oxid Med Cell Longev. 2022 Apr 19;2022:4608914.

See more customer validations on www.MedChemExpress.com

### **REFERENCES**

- [1]. Landry DW, et al. Epithelial chloride channel. Development of inhibitory ligands. J Gen Physiol. 1987 Dec;90(6):779-98.
- [2]. Landry DW, et al. Purification and reconstitution of chloride channels from kidney and trachea. Science. 1989 Jun 23;244(4911):1469-72.
- [3]. Xiaoqin Fan, et al. A homogeneous time-resolved fluorescence-based high-throughput screening for discovery of inhibitors of Nef-sdAb19 interaction. Int J Oncol. 2015 Oct;47(4):1485-93.

 $\label{lem:caution:Product} \textbf{Caution: Product has not been fully validated for medical applications. For research use only.}$ 

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA