


(24S)-24,25-Dihydroxyvitamin D3

Cat. No.:	HY-15439
CAS No.:	55700-58-8
Molecular Formula:	C ₂₇ H ₄₄ O ₃
Molecular Weight:	416.64
Target:	VD/VDR
Pathway:	Vitamin D Related/Nuclear Receptor
Storage:	-20°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)



SOLVENT & SOLUBILITY

In Vitro	DMSO : 100 mg/mL (240.02 mM; Need ultrasonic)				
		Solvent Concentration	Mass		
	Preparing Stock Solutions			1 mg	5 mg
		1 mM		2.4002 mL	12.0008 mL
		5 mM		0.4800 mL	2.4002 mL
	10 mM		0.2400 mL	1.2001 mL	
Please refer to the solubility information to select the appropriate solvent.					
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (6.00 mM); Clear solution				
	2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: 2.5 mg/mL (6.00 mM); Suspended solution; Need ultrasonic				
	3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (6.00 mM); Clear solution				

BIOLOGICAL ACTIVITY

Description	(24S)-24,25-Dihydroxyvitamin D3 ((24S)-24,25-Dihydroxycholecalciferol) is an inactive form of vitamin D3 which undergoes various levels of hydroxylation to form active vitamin D3 analogs.
In Vitro	(24S)-24,25-Dihydroxyvitamin D3 ((24S)-24,25-Dihydroxycholecalciferol) is a isomer of 24,25-dihydroxyvitamin D3. Vitamin D is metabolized in the liver to 25-hydroxyvitamin D3, and in the kidney to 1,25-dihydroxyvitamin D3 and 24,25-dihydroxyvitamin D3 ^[5] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

- [1]. Yamada S, Ino E, Takayama H et al. Differences in the side-chain metabolism of vitamin D3 between chickens and rats. *Proc Natl Acad Sci U S A*. 1985 Nov;82(22):7485-9.
- [2]. Kobayashi N, Higashi T, Saito K et al. Specificity of polyclonal antibodies raised against a novel 24,25-dihydroxyvitamin D3-bovine serum albumin conjugant linked through the C-11alpha or C-3 position. *J Steroid Biochem Mol Biol*. 1997 May;62(1):79-87.
- [3]. Larsson D, Nemere I, Sundell K. Putative basal lateral membrane receptors for 24,25-dihydroxyvitamin D(3) in carp and Atlantic cod enterocytes: characterization of binding and effects on intracellular calcium regulation. *J Cell Biochem*. 2001 Aug 1-9;83(2)
- [4]. Nemere I, Yazzie-Atkinson D, Johns DO, Larsson D. Biochemical characterization and purification of a binding protein for 24,25-dihydroxyvitamin D3 from chick intestine. *J Endocrinol*. 2002 Jan;172(1):211-9.
-

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA