Proteins

Product Data Sheet

Ertapenem

Cat. No.: HY-A0294 CAS No.: 153832-46-3 Molecular Formula: $C_{22}H_{25}N_3O_7S$

Molecular Weight: 475.51

Target: Antibiotic; Bacterial Pathway: Anti-infection

Please store the product under the recommended conditions in the Certificate of Storage:

BIOLOGICAL ACTIVITY

Description Ertapenem (MK-0826) is a broad spectrum and long acting β-lactam antibiotic. Ertapenem has a broad-spectrum antianaerobic activity against a variety of anaerobes with a mode MIC of 0.12 µg/mL. Ertapenem can be used for the research of severe infections caused by bacteria in the skin, lungs, stomach, pelvis, and urinary tract^{[1][2]}.

IC₅₀ & Target **B-lactam**

In Vitro

Ertapenem (0-100 $\mu g/mL$ approximately, 48 h) is active against 99.1% of all anaerobes with a mode MIC of 0.12 $\mu g/mL$ and MIC_{90} of 1 µg/mL, and MIC's ≥ 8 µg/mL for B.fragilis and B.vulgatus species, respectively^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Cell Viability Assay^[1]

Cell Line:	B. fragilis (ATCC 25285), B. thetaiotaomicron (ATCC 29741), and Eubacterium lentum (ATCC 43055)
Concentration:	0-100 μg/mL approximately
Incubation Time:	48 h
Result:	Inhibited 99.1% of all isolate with a mode MIC of 0.12 μ g/mL and MIC ₉₀ of 1 μ g/mL, and 98.8% of the isolates were susceptible among the B. fragilis group.

In Vivo

Ertapenem (Subcutaneous injection, 0-10 mg/kg, 0-120 h after infection, S. aureus thigh tissue infection model) shows > 3 log₁₀ CFU reduction of organism at 10 mg/kg, and maintains the activity with 3.3 and 4.4 log₁₀ CFU eliminated at 2 mg/kg^[2]. Ertapenem (Subcutaneous injection, 4h after infection, systemic infection model) is active against all gram-positive organisms, and is also active against gram-negative organisms tested with ED₅₀s of <0.25 mg/kg/dose^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Model:	S. aureus thigh tissue infection model (DBA/2 mice) $^{[2]}$
Dosage:	0.5,1, 2, 5, 10 mg/kg (given at 2, 6, 10, 24, 48, 72, 96, 120 h)
Administration:	Subcutaneous injection (0.5 mL after infection)

Result:	Displayed > 3 \log_{10} CFU reduction of organism compared to non-antibiotic-treated controls at 10 mg/kg.
	Maintained the activity with 3.3 and 4.4 \log_{10} CFU eliminated at 2 mg/kg.
Animal Model:	Systemic infection model (DBA/2 female mice, viral antibody-free CD-1 female mice) ^[2]
Dosage:	0-3 mg/kg approximately
Administration:	Subcutaneous injection (0.5 mL, begin immediately and 4 h after infection)
Result:	Showed activity against all gram-positive organisms, and also ram-negative organisms tested with ED $_{50}$ s of <0.25 mg/kg/dose.
Animal Model:	CD-1 mice, rats ^[2]
	10 mg/kg approximately
Dosage: Administration:	Intraperitoneal injection (pharmacokinetic assay)
Result:	Exhibited an AUC _{0-∞} ranging from 1.8-21.82 μg•hr/mL in tissue in mice following a 10-mg/kg i.p. dose.
	Exhibited slow clearance rate with a $t_{1/2\beta}$ of 3.2 h, Clp of 0.47 mL/min/kg, AUC ₀₋₈ of 284.15 μ g•hr/mL.

CUSTOMER VALIDATION

- Nat Commun. 2022 Mar 2;13(1):1116.
- Nat Commun. 2021 Jul 22;12(1):4461.
- J Antimicrob Chemother. 2023 Jul 31;dkad229.
- J Antimicrob Chemother. 2020 Jul 1;75(7):1850-1858.
- Biomed Res Int. 2018 Jul 2;2018:3579832.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Kenneth E Aldridge. Ertapenem (MK-0826), a new carbapenem: comparative in vitro activity against clinically significant anaerobes. Diagn Microbiol Infect Dis. 2002 Oct;44(2):181-6.

[2]. C J Gill, et al. In vivo activity and pharmacokinetic evaluation of a novel long-acting carbapenem antibiotic, MK-826 (L-749,345). Antimicrob Agents Chemother. 1998 Aug;42(8):1996-2001.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898 Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Page 2 of 2 www.MedChemExpress.com